


TECHNICAL DRAWING



Art.	Vaso	A	B	C	D	E	Apertura frontale	Coperchio
440	470; 470-V	525	385	360	260	155	SI	SI

Misure in mm.

PHOTO



**PRODUCT FAMILY:** SANITARYWARE

**CODE:** 440

**DESCRIPTION:** Toilet seat with lid for toilet artt. 470, 470-B, 470V, 470V-B.

**APPROVALS/CERTIFICATIONS:**

TUV system certification UNI EN ISO 9001:2008.

**PRODUCT FEATURES:**

Toilet seat with lid for toilet artt. 470, 470-B, 470V, 470V-B.. Chipboard core coated with cast polyester (thick. 2mm). Polyester coating, compared to superficial painting, offers a better resistance to wear, therefore an increased duration of the article. Complete with chrome-plated brass hinges. Expansion fixing.

NB: EACH SEAT FITS A SPECIFIC MODEL OF TOILETS. IF MOUNTED ON OTHER TOILETS, IT COULD NOT ADAPT PERFECTLY, THUS

CAUSING THE HINGES TO BREAK. IT IS HIGHLY RECOMMENDED THAT YOU USE THE SEAT CORRESPONDING TO EACH TOILET. STANDARD SEATS ARE USED FOR TOILETS WITH MORE COMMON SHAPES. NEVERTHELESS, WE RECOMMEND CHECKING THE SHAPE AND DIMENSIONS BEFORE THE INSTALLATION.

**FIXING/INSTALLATION:**

Fasten the seat to the toilet using the included screws, making sure that the fasteners are securely tightened to avoid backlashes. The front bumpers must lean against the toilet edge, the back ones must be 3-4 mm raised.

**MAINTENANCE:**

Check the correct fastening of the seat periodically. Avoid backlashes to the fasteners and incorrect movements, which could damage the hinges. Clean with water and soap or other suitable products, do not use abrasive or corrosive agents. Dry with a soft cloth.

**PRODUCT COMPOSITION:**

N° 1 Seat with lid - Chipboard, Cast Polyester  
N° 2 Hinges - Chrome-plated brass  
N° 6 Bumpers - Linear polyethylene  
N° 2 Lid bumpers - Linear polyethylene

**PACKAGING OF THE PRODUCT:**

Packed in cardboard box 380 x 555 x 55 mm.  
Multiple packing (5 pcs) in individual cardboard boxes + 580 x 330 x 400 mm box.  
Weight: 3,7 kg

The above mentioned information cannot be considered as binding. Thermomat Saniline Srl reserves the right to make technical and aesthetic changes without any prior notice.

ST_PRODOTTO: 440 - 1	APPROVATO R.T.
REV.00 DEL 25/09/09	

**THERMOMAT SANILINE SRL**  
Via M. Curie, 7 - 46034 Borgo Virgilio (Mantova) ITALY  
Tel. + 39 0376 4473  
www.thermomat.com - e-mail: info@thermomat.com

